



Scientific / Metrology Instruments  
Energy-dispersive X-ray fluorescence spectrometer

Solutions for Innovation

# JSX-1000S

# ElementEye

The Smart Solution for Elemental Analysis



JEOL Ltd.

## Smarter solutions for more users

The JSX-1000S is an X-ray fluorescence spectrometer that provides easy, quick elemental analysis using touch panel operation.



# High performance with simple operation

Based on many years of experience, JEOL has incorporated proprietary algorithms in the Fundamental Parameter (FP) method. This FP method delivers highly accurate qualitative and quantitative analysis results, without using a standard sample.

JEOL's own SDD (silicon drift detector) and newly-developed optical system, in combination with filters designed to handle the entire energy range, make it possible to achieve high-sensitivity analysis.

The sample chamber vacuum unit (option) further increases detection sensitivity for lighter elements.

## Simple Operation

*Just Set the sample and Touch the screen.*

Analysis values and pass/fail assessment results are updated in real-time, and the spectrum can be checked easily.



Swipe the screen to switch the display, and check the spectrum


## Sensitive analysis throughout the entire energy range

High-sensitivity analysis can be performed across the entire energy range using a maximum of 9 types of filters and a sample chamber vacuum unit.



\*Cl, Cu, Mo, and Sb filters are optional


Table 1 Element correspondence with each filter and sample chamber vacuum unit.



## Options

### • Sample chamber vacuum unit

Performing analysis in a vacuum makes it possible to improve sensitivity for lighter elements (Na ~ Al).



### • Auto sample changer

A maximum of 12 samples can be exchanged automatically, allowing automatic continuous acquisition.



## High-precision quantitative analysis

The new Smart FP (Fundamental Parameter) method makes it possible to obtain highly-accurate quantification results without the need to prepare a standard sample, and includes automatic correction for thickness and residual ingredient balance.

(The residual balance correction and thickness correction functions are only applicable to organic samples)


Table 2 Example : Polyethylene

Thickness	Correction	Cr	Zn	Cd	Pb	Automatic balance
0.5mm	No	0.008	0.037	0.001	0.002	99.76
3.8mm		0.012	0.109	0.004	0.006	99.64
0.5mm	Yes	0.011	0.137	0.015	0.010	99.54
3.8mm		0.011	0.134	0.016	0.011	99.55
Standard value		0.010	0.125	0.014	0.010	

## High Throughput

The new high-sensitivity SDD (silicon drift detector) and short-path optical system achieve throughputs up to 7.6 times faster than the previous JEOL system (JSX-3100R II)

Table 3 Comparison of Cd 10 ppm / Pb 50 ppm measurement times



## Providing Solutions

With solution based applications, the desired analysis can be executed automatically according to pre-recorded recipes. Simply select the desired solution icon from the solution application list for automated analysis and display of results.

..... *Solution applications offer simplified analysis for a wide range of industries* .....

### Standard Solution

- RoHs analysis (Metals)
- Simple analysis (Metals: Air / Vacuum\*)
- Simple analysis (Oxides: Air / Vacuum\*)
- RoHs analysis (Metals)
- Simple analysis (Organic materials: Air / Vacuum\*) (\*option)

### Optional Solution

- Ni plating screening
- Halogen screening
- Sn plating screening



# Smart Assist from ElementEye




With an X-ray fluorescence spectrometer you can perform elemental analysis with little or no sample preparation. The easy operation of ElementEye assists your analysis work.

## SEM-EDS

XRF solid sample analysis assist



SEM: Backscattered electron image





Zirconium (Zr) is found in sand from ocean and rivers; but not all sand contains Zr. Using XRF makes it possible to quickly determine whether there is any Zr in a sand sample. Using SEM-EDS allows you to identify which grains of sand contain Zr. Combining XRF and SEM-EDS makes it easy to pinpoint the location of a target element.

## XPS XRF thin-film analysis assist



Printed circuit board layers




1	1, Cu
2	2, Ni
3	3, Cr
4	4, Polyimide



When acquiring a depth profile of a multi-layer sample using XPS, knowing the thickness of each layer in advance makes it possible to set the optimal etching conditions. With XRF it is easy to quickly determine the thickness of each layer. The depth profile data on the right was obtained by setting the etching conditions based on the layer thicknesses determined with XRF. With XPS, in addition to the film structure of the main components, it is also possible to identify non-uniform distributions of oxygen.

## NMR

XRF fluid sample analysis assist



NMR is widely used for the structural analysis of organic materials, measuring mainly hydrogen ( $^1\text{H}$ ) and carbon ( $^{13}\text{C}$ ). For measurements of nuclei other than  $^1\text{H}$  and  $^{13}\text{C}$ , the sensitivity varies depending on the concentration of the species. Therefore, obtaining composition information for the sample beforehand enables estimation of the appropriate acquisition conditions. Using XRF it is easy to determine the types and concentrations of elements, like Pt, contained in a sample.

# Main Specifications & Options

## Specifications

Detectable element range	Mg ~ U Na ~ U (option)
X-ray generator	5~50 kV 1mA
Target	Rh
Filters: Maximum of 9 types, Automatic exchange	OPEN, ND ((neutral density filter), Cr, Pb, Cd Cl, Cu, Mo, Sb (options)
Collimators: 3 types, Automatic exchange	0.9mm, 2mm, 9mm
Detector	Silicon drift detector
Sample chamber size	300mm $\phi$ x 80mmH
Sample chamber atmosphere	Air atmosphere / Vacuum (option)
Sample chamber observation mechanism	Color camera
Computer for operation	Windows <sup>®</sup> desktop computer with touch panel
Analysis software	Qualitative analysis (Auto qualitative, KLM marker, Sum peak display, Spectrum search) Quantitative analysis (Bulk FP method, RoHS solution (Cd, Pb, Cr, Br, Hg) Simple analysis solution Report generation tool
Daily check software	Tube aging, Energy check, Intensity check

\* Windows<sup>®</sup> is a registered trademark of Microsoft<sup>®</sup> Corp. in the U.S.A. and other countries

## Main options

- Sample chamber vacuum unit
- Auto sample changer
- Filter set
- Filter FP method analysis software
- Thin-film FP method analysis software
- Sum peak removal software
- Ni plating screening solution
- Sn plating screening solution
- Halogen screening solution

## Installation requirements


Power supply	Single phase 100 to 120 V $\pm$ 10%, 15 A, 50 / 60 Hz, Type D grounding
Environment	Temperature: 18 to 28 $^{\circ}$ C Humidity: 30 to 70% (no condensation) Not subject to corrosive gas, dust, or discernible vibration Flows of heated or cooled air should not be directed at the system

\*Leave a maintenance space of at least 50 cm on all sides of the unit

\*The external appearance and specification are subject to change without notice

\*The customer must provide the work table and chair used for this instrument and PC.

## JSX-1000S installation layout



\*Legal documentations are required for instruments using X-rays, according to the regulations of countries where you are using them.

## World Wide Service and Support

At JEOL, our goal is customer satisfaction. To meet customer requirements, JEOL has made great efforts to produce a high performance instrument with simple operation in the Element Eye.

JEOL has 20 subsidiary companies and a network of service support centers all over the world to maximize the up-time of your instrument.

Furthermore, our application engineers are available to provide you with technical information ranging from basic elemental analysis to technical solutions required for operation and maintenance, at local service centers around the world.

JEOL will continue to provide you with technical information, qualified repairs, consumable parts and attachments on a long-term basis.



## World Wide Support Network



High technology for quality assurance and the environment.

Certain products in this brochure are controlled under the "Foreign Exchange and Foreign Trade Law" of Japan in compliance with international security export control. JEOL Ltd. must provide the Japanese Government with "End-user's Statement of Assurance" and "End-use Certificate" in order to obtain the export license needed for export from Japan. If the product to be exported is in this category, the end user will be asked to fill in these certificate forms.

\*specifications subject change without notice.

**JEOL** JEOL Ltd.

1-2 Musashino 3-chome Akishima Tokyo 196-8558 Japan Sales Division Telephone: +81-3-6262-3560 Facsimile: +81-3-6262-3577

<http://www.jeol.com/>

•AUSTRALIA & NEW ZEALAND/ JEOL(AUSTRALASIA) Pty.Ltd. Suite 1, L2 18 Aquatic Drive - Frenchs Forest NSW 2086 Australia •BELGIUM/ JEOL (EUROPE) B.V. Planet II, Gebouw B Leuvensesteenweg 542, B-1930 Zaventem Belgium •CANADA/ JEOL CANADA, INC. 3275 1ere Rue, Local #8 St-Hubert, QC J3Y-8Y6, Canada •CHINA/ JEOL(BEIJING) CO., LTD. Room B1010/1110, Wantong New World Plaza No. 2 Fuchengmenwai Street, Xicheng District, Beijing 100037, P.R.China •EGYPT/ JEOL SERVICE BUREAU 3rd Fl. Nile Center Bldg., Nawal Street, Dokki, (Cairo), Egypt •FRANCE/ JEOL (EUROPE) SAS Espace Claude Monet, 1 Allée de Giverny 78290, Croissy-sur-Seine, France •GERMANY/ JEOL (GERMANY) GmbH Oskar-Von-Miller-Strasse 1a, 85386 Eching, Germany •GREAT BRITAIN & IRELAND/ JEOL (U.K.) LTD. JEOL House, Silver Court, Watchmead, Welwyn Garden City, Herts AL7 1LT, U.K. •INDIA/ JEOL INDIA Pvt. Ltd. Elegance Tower, Level 2, 212B Old Mathura Road, Jasola Business District Centre, Near Apollo Hospital Jasola, New Delhi 110 025, India •ITALY/ JEOL (ITALIA) S.p.A. Palazzo Pacinotti - Milano 3 City, Via Ludovico il Moro, 6/A 20080 Basiglio(MI) Italy •KOREA/ JEOL KOREA LTD. Dongwoo Bldg. 7F, 1443, Yangjae Daero, Gangdong-Gu, Seoul, 134-010, Korea •MALAYSIA/ JEOL(MALAYSIA) SDN.BHD. 508, Block A, Level 5, Kelana Business Center, 97, Jalan SS 7/2, Kelana Jaya, 47301 Petaling Jaya, Selangor, Malaysia •MEXICO/ JEOL DE MEXICO S.A. DE C.V. Arkansas 11 Piso 2 Colonia Napoles Delegacion Benito Juarez, C.P. 03810 Mexico D.F., Mexico •SCANDINAVIA/ SWEDEN JEOL (Skandinaviska)AB: Hammarbacken 6A, Box 716, 191 27 Sollentuna Sweden •SINGAPORE/ 2 Corporation Road #01-12 Corporation Place Singapore 618494 •TAIWAN/ JIE DONG CO., LTD. 7F, 112, Chung Hsiao East Road, Section 1, Taipei, Taiwan 10023 Republic of China •THE NETHERLANDS/ JEOL (EUROPE) B.V. Lireweg 4, NL-2153 PH Nieuw-Vennep, The Netherlands •USA/ JEOL USA, INC. 11 Dearborn Road, Peabody, MA 01960, U.S.A.

•Please confirm other territories by the web site.

No.1301H445C Printed in Japan, Bn