



ratiopetta®

Pipetting is one of the most frequently used applications in laboratories. This is not only where precision and reliability are important, users also need products that are ergonomic, light and convenient. The new ratiopetta® from ratiolab® meets these demands while offering outstanding value for money.

Product properties

- obtainable in volumes from 0.1 µl – 5000 µl
- full autoclavable, 20 min, 121°C
- universal cone for standard pipet tips
- easy volume setting by smooth-running setting ring
- fatigue-free work possible for the analyst over longer pipetting series
- touch-sensitive pressure points for pipetting and discharging
- modern, innovative materials
- individual checking according to ISO8655/DIN12650
- delivery includes a test certificate
- individual serial numbers
- suitable for right- and left-handers
- easy to calibrate
- corrosion-resistant piston and tip ejector
- 1 year guarantee

Delivery includes:

- pipet
- certificate of calibration
- operating instructions
- bench rack
- calibration key
- lubricant

Matching pipet tips can be found on pages 20-25.





Single-channel Precision Pipet ratiopetta®

CE IVD

Volume Range µl	Accuracy* < ± %	Precision* < %	Quantity per Pack	Order No.
0.1 – 2.5	3.00	3,00	1	47 80 025
0.5 – 10	1.50	1.50	1	47 80 010
2.0 – 20	1.20	1,00	1	47 80 020
10 – 100	1.00	0.40	1	47 80 100
20 – 200	0.80	0.40	1	47 80 200
100 – 1000	0,70	0.25	1	47 81 000
1000 – 5000	0.60	0.30	1	47 85 000

8-channel Precision Pipet ratiopetta®

CE IVD

Volume Range µl	Accuracy* < ± %	Precision* < %	Quantity per Pack	Order No.
0.5 – 10	2.50	2.50	1	47 80 018
5.0 – 50	1.50	1.00	1	47 80 058
50 – 300	1.00	0.50	1	47 80 308

12-channel Precision Pipet ratiopetta®

CE IVD

Volume Range µl	Accuracy* < ± %	Precision* < %	Quantity per Pack	Order No.
0.5 – 10	2.50	2.50	1	47 80 012
5.0 – 50	1.50	1.00	1	47 90 052
50 – 300	1.00	0.50	1	47 90 302

Matching pipet tips can be found on pages 20-25.

Pipettor stands ratiopetta®

- for the clearly laid out storing of up to 5 pipettors of different brands with finger rest
- compatible with most multi-channel-pipettors with finger rest

Product	Quantity per Pack	Order No.
pipettor stand for ratiopetta®	1	94 20 024

* Three measurements were carried out for all the pipets: minimum, maximum and mean volumes. The data shown are based on the mean volumes. On request, we can provide you with all of the results that were measured for your pipet. Please let us know the order number and the serial number of your Ratiopetta.