









Dispenser Tips

	Volume setting				
	1	2	3	4	5
 1	1	2	3	4	5µl
 2	10	20	30	40	50µl
 3	25	50	75	100	125µl
 4	50	100	150	200	250µl
 5	0.25	0.5	0.75	1.0	1.25ml
 6	0.5	1.0	1.5	2.0	2.5ml
 7	1	2	3	4	5ml
 8					

Dispenser Tips

CE IVD

For serial pipettings

Dispenser Tips fit onto ratiolab® dispettea and Eppendorf Multipipette® 4780, Ritter Dispenser and Distriman (Gilson). Not autoclavable.

8 Adapters for 25 ml and 50 ml Dispenser Tips must be ordered separately.

Volume ml	Fig. No.	Description	Quantity per Pack	Order No.
0.5	1	non-sterile	100	29 00 101
1.25	2	non-sterile	100	29 00 102
2.5	3	non-sterile	100	29 00 103
5.0	4	non-sterile	100	29 00 104
12.5	5	non-sterile	100	29 00 105
25	6	non-sterile	50	29 00 109
50	7	non-sterile	25	29 00 106
Adapter	8	non-sterile	10	29 00 107

0.5	1	sterile, individually packed	100	29 10 101
1.25	2	sterile, individually packed	100	29 10 102
2.5	3	sterile, individually packed	100	29 10 103
5.0	4	sterile, individually packed	100	29 10 104
12.5	5	sterile, individually packed	100	29 10 105
25	6	sterile, individually packed	25	29 10 109
50	7	sterile, individually packed	25	29 10 106
Adapter	8	sterile, individually packed	5	29 10 107

Dispetta Dispenser



Dispetta Pipettor for Serial Dispensing

CE IVD

Convenient and robust instrument for serial dispensing with Dispenser Tips.

Five volume settings combined with seven Dispenser Tips sizes result in 35 different volumes from 10 µl to 5 ml.

The Dispenser Tip is filled similar to a syringe. After preselecting the dispensing volume the pipet button has to be pressed for each dispensing step.

- ergonomic design and light weight for convenient handling
- smooth and fatigue-free operation even with long series
- precise and convenient volume setting
- clearly visible volume adjustment, optimized readability
- table with volume scale fixed directly onto the instrument
- precise gear mechanism for accurate and reproducible pipettings
- safe, rapid and easy exchange of Dispenser Tips
- maintenance-free
- long life due to robust design
- manual dispenser without electronic parts
- 2 years guarantee

Product	Quantity per Pack	Order No.
ratiolab® dispetta dispenser	1	29 80 200